
LIVO

Item number ottobock. LIVO: 476C00 = SK*; 477C00 = SV*

DE	Bedienungsanleitung	2
GB	Instructions for use	3
NL	Gebruiksaanwijzing	4
CZ	Návod k obsluze	5
ES	Instrucciones de uso	6
IT	Istruzioni per l'uso	7
FR	Mode d'emploi	8
GR	Οδηγίες χειρισμού	9
PL	Instrukcja obsługi	10
PT	Instruções de operação	11
HR	Upute za rukovanje	12
NO	Bruksanvisning	13
FI	Käyttöohjeiden	14
DK	Betjeningsvejledning	15
RU	Инструкция по эксплуатации	16
SE	Bruksanvisningar	17
SI	Navodila za uporabo	18
HU	Üzemeltetési utasítások	19
SK	Návod na obsluhu	20
KR	사용 설명서	21
CN	操作说明	22
CN	操作說明 (Traditional Chinese)	23

Dekubitusprophylaxe Sitzkissen

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Decubitus prophylaxis seat cushion

Please read the instructions and the notes on safety carefully before using the system to ensure safe operation!

Product name	ottobock. LIVO
Scope of delivery	decubitus prophylaxis seat cushion with cover and Instructions of use
Indication	Decubitus prophylaxis
Contraindication	None
Range of application	Decubitus prophylaxis
Usage	The seat cushion is used in wheelchairs. The cushion size should be selected according to the individual wheelchair dimensions.
Patient weight min. – max.	0 – 125 kg
Material	PU-foams (HR/Visco); PU-coating/ /Mesh cover
Sizes	See product label
Weight of product	Approx. 0,5 kg
Medical device class 1 according 93/42 EWG and medical product law (MPG)	
Technical safety / Flame retardancy	Flammability test according to DIN EN 1021-1 & EN 1021-2
Biocompatibility	DIN EN ISO 10993-5, DIN EN ISO 10993-1
Notes on safety	Improper use Risk of damage to the cushion surface due to user error • Do not allow the wheelchair cushion to come into contact with sharp objects. This applies to animals such as pet cats with sharp claws. • Do not expose the wheelchair cushion to excessive heat or open flame. • If the wheelchair cushion is expected to come into contact with liquid, such as spilt drink episodes of incontinence, always use it in conjunction with a liquid-repellent cover. Contact your therapist or specialist dealer to obtain a spare Ottobock cover.
Guaranty	3 years - Guarantee applies to functional properties of the cushion core (pro rata guarantee).
Washing instructions (cover)	

Cleaning & Maintenance Methods

Home use	Clinical use
While patient using the product	
Cover and foam core can be wiped with commercial disinfectants made by isopropanol or other DHGM listed products; or Chemothermal cleaning of the cover.	Cover can be cleaned by chemothermal cleaning or with commercial disinfectants made by isopropanol or other DHGM listed products.
On every change of patient / Product repeated using	
Cover can be wiped with commercial disinfectants made by isopropanol or other DHGM listed products; or chemothermal cleaning of the cover.	Cover can be cleaned by chemothermal cleaning or with commercial disinfectants made by isopropanol or other DHGM listed products.
Autoclaving the foam core according to RKI-guidelines; or the foam core without cover can be cleaned with household detergents. This can be achieved by saturating the contaminated site with a disinfectant (e.g. Ecolab Oxyfoam, consisting of hydrogen peroxide). The unfilled foam core may be sanitized with a common household disinfectant (e.g., Bacillol) by soaking the contaminated site or spraying the entire foam core. Then let the foam core dry in the air.	Autoclaving the foam core according to RKI-guidelines.
No further maintenance necessary. This product can repeated using.	

Instructions for use

1. Take the cushion out of the box.
2. Put the cushion onto the wheelchair or seat. The ottobock label must be found on the front right.

RKI = Robert-Koch-Institut (RKI) is a part of the Federal Ministry of health in Germany. It is responsible for disease prevention and control.
DGHM = deutsche Gesellschaft für Hygiene und Mikrobiologie (German Society for Hygiene and Mikrobiology)
MPG = Medizinproduktegesetz (Medical Devices Act)
Autoclaving = fractionated vacuum method at 105 degrees C and 55 minutes

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Otto Bock Mobility Solutions GmbH
Funke Medical AG